




# Disposables and Accessories Catalogue



[www.laborie.com](http://www.laborie.com)

## CATALOGUE INDEX

PAGE	PART #	DESCRIPTION	UNITS/BOX
<b>LABORIE CATHETERS</b>			
1	CAT206	Twin lumen catheter for cystometry 6 FR	10
1	CAT208	Twin lumen catheter for cystometry 8 FR	10
1	CAT307	Triple lumen catheter for UPP 7 FR	10
1	CAT411	Rectal balloon catheter 10 FR (latex)	10
2	CAT506	Pediatric Rectal balloon catheter 8 FR (silicone)	10
2	CAT509	Rectal balloon catheter 10 FR (silicone)	10
<b>AIR CHARGED CATHETERS</b>			
3	CAT875	TDOC, Abdominal catheter, 7 FR	10
3	CAT895	TDOC, Single Sensor catheter, 7 FR	10
3	CAT880	TDOC, Dual Sensor catheter, 7 FR	10
3	DIS500	Stylet 10	1
4	DIS950	Velcro Stabilizer Leg Strap	1
4	DIS450	Air Charged Calibration Chamber	1
<b>UDS TUBING</b>			
5	TUB101	Pressure tubing set with stopcock	25
5	TUB500	Laborie pump tubing set	25
5	TUB501	Laborie pump head tubing	25
5	TUB219	Perfusion line with 2 line connections	5
<b>UDS TRANSDUCER CARTRIDGES AND DOMES</b>			
6	DIS102	Luer Lock Plug	100
6	DIS103	Nova Dome male/female (UDS 64 based systems)	50
6	DIS101	Nova Dome male/male (UDS 64 based systems)	50
6	DIS130	Pressure cartridge with luer lock plug (UDS94 based systems)	25
<b>EMG</b>			
7	ELE425	EMG patch electrodes, 3/pk, without lead cable	100
7	ELE428	EMG patch electrodes, 3/pk, with lead cable	100
7	ELE350	EMG lead set non shielded: 2 red & 1 green	1
7	ELE200	EMG lead set shielded: 2 red & 1 green	1
7	CAB154	EMG Cable, Non Shielded	1
7	CAB155	EMG Cable, Shielded	1
<b>BIOFEEDBACK/STIMULATION</b>			
8	PRB450	Vaginal Classic Probe	1
8	PRB055	Anal Short Probe	1
8	PRB085	Vaginal Atrophic Probe	1
8	PRB425	Vaginal Tampon Probe	1
9	PRB113	Vaginal Manometry Probe	1
9	PRB108	Anal Manometry Probe	1
9	PRB080	Tubing for Vaginal Manometry Probe	10
9	PRB078	Tubing for Anal Manometry Probe	10
10	ELE625	Patch Electrodes (4pk)	15
<b>MISCELLANEOUS</b>			
11	DIS173	Graduated beaker	10
11	DIS155	Thermal printer paper rolls	5
11	DIS302	Pressure cuff	1
11	CHA171	Standard Commode Chair	1
12	CHA181	Folding Commode Chair	1
12	CHA172	Short Funnel	1
12	CHA876	Long Funnel	1
12	STA700	Adjustable Uroflow Stand	1

### GENERAL SALES CONDITIONS:

#### Only complete boxes are sold

Terms of delivery: 14 days to 6 weeks depending on the quantity ordered.

Items are Ex-works and excl. VAT or other applicable custom taxes or charges.


Please specify the custom broker if applicable.

Written order requested by email or fax. Order number or other reference is required.

## Laborie Twin Lumen, 6Fr

---


Part Number: CAT206



## Laborie Twin Lumen, 8Fr

---


Part Number: CAT208



## Laborie Triple Lumen, 7Fr

---


Part Number: CAT307



## Laborie Rectal Balloon, 10Fr (latex)


---

Part Number: CAT411




## Laborie Rectal Balloon Pediatric, 8Fr (silicone)

Part Number: CAT506



## Laborie Rectal Balloon, 10Fr (silicone)

Part Number: CAT509



## Abdominal Catheter (rectal or vaginal), 7Fr

Single sensor, no lumen


**Part Number:** CAT875



## Bladder Catheter, 7Fr

Single sensor, with lumen


**Part Number:** CAT895



## Bladder/Urethral Catheter, 7Fr

Dual sensor, with lumen


**Part Number:** CAT880



## Stylet 10

For use with CAT895

**Part Number:** DIS500



\*Urodynamic system must be equipped specifically for use with TDOC catheters

## Velcro Leg Strap

---

Part Number: DIS950



## Calibration Chamber

---

Part Number: DIS450



\*Urodynamic system must be equipped specifically for use with TDOC catheters

## Laborie Pressure Tubing with Stopcock

For connecting the transducer to the catheter

**Part Number:** TUB101



## Laborie Pump Tubing Set

**CAUTION:** This tubing set must be used with all Laborie Urodynamics systems to assure accuracy of the infused volume to the bladder.

**Part Number:** TUB500



## Laborie Pump Head Tubing

Laborie pump head tubing

**Part Number:** TUB501



## Perfusion Tubing

Perfusion line with 2 line connections

**Part Number:** TUB219

**Note:** Photo of product not shown

## Luer Lock Plug

---

**Part Number:** DIS102



## Novadome Male/Female

---

For use with UDS64 systems

**Part Number:** DIS103

**NOTE:** For Novadome with male/male connectors order  
**Part Number:** DIS101 (not shown here)



## Pressure Transducer Cartridge with Luer Plug

---

For use with all current UDS94 systems

**Part Number:** DIS130



### EMG Gel Patch Electrodes

---

Snap pre-gelled EMG Electrodes, 3 per package

**Part Number:** ELE425




### EMG Gel Patch Electrodes w/lead cables

---

EMG Electrodes with lead cables, 3 per package

**Part Number:** ELE428




### EMG Leads - Non Shielded

---

**Part Number:** ELE350

**NOTE:** for video UDS systems order the shielded EMG leads,  
**Part Number:** ELE200 (not shown here)




### EMG Cable, Non Shielded

---

For use with all non video UDS systems

**Part Number:** CAB154

**NOTE:** for video UDS systems order the shielded EMG cable,  
**Part Number:** CAB155 (not shown here)



**Vaginal Classic Probe**

---

Part Number: PRB450



**Anal Short Probe**

---

Part Number: PRB055



**Vaginal Atrophic Probe**

---

Part Number: PRB085



**Vaginal Tampon Probe**

---

Part Number: PRB425



**Vaginal Manometry Probe**

---

Part Number: PRB113



**Anal Manometry Probe**

---

Part Number: PRB108



**Tubing For PRB113 Vaginal Manometry Probe**

---

Part Number: PRB080



**Tubing For PRB108 Anal Manometry Probe**

---


Part Number: PRB078



**Surface Patch Electrodes, 4 per package**

---

**Part Number:** ELE625



## Graduated Beaker

---

Part Number: DIS173



## Thermal Printer Paper

---


Part Number: DIS155



## Pressure Cuff

---

Part Number: DIS302



## Standard Commode Chair

---

Part Number: CHA171



## Folding Commode Chair

---

**Part Number:** CHA181



## Short Funnel

---

For use with all Urocap and all UDS64 based systems

**Part Number:** CHA172



## Long Funnel

---

For use with all current UDS94 systems

**Part Number:** CHA876



## Adjustable Uroflow Stand

---

Funnel not included

**Part Number:** STA700



## **AQUARIUS TT**™

The ultimate in UDS technology utilizing our sleek patient tower with hi-speed EMG sampling up to 5000 samples per second, electro-mechanical height controlled transducers and a powerful processor. Aquarius TT includes the most comprehensive software package and a custom built industrial computer. Perform Video Urodynamics studies with our state-of-the-art Smart Capture™ technology. The Aquarius TT is the right choice for the most sophisticated urologist or research institution.



## **TRITON**™

Our state-of-the-art wireless UDS system with a patient tower that can be located close to the patient while the PC, tablet or laptop can be located anywhere in the room. Hi-speed EMG sampling up to 3000 samples per second, electro-mechanical height controlled transducers and a powerful processor. A broad range of software features are standard with additional options available including an upgrade to perform video urodynamics.



## **DORADO KT**™

The latest innovation in UDS systems by Laborie. All components are contained in an attractive, ergonomically-sound mobile cart. Optional software packages can be customized to suit your needs. Adjustable shelves and workspace trays allow for customizable configurations. The KT cart provides a high degree of in-office portability with lockable wheelbase.



## **DELPHIS™**

Laborie led the industry in 2002 introducing the Delphis utilizing wireless Bluetooth™ technology. Compact and affordable, it is an excellent system for a clinic. The Delphis has a wireless connection between the laptop and UDS processor which allows the user to remove the laptop from the cart during patient voiding for privacy while viewing the live study.



## **DELPHIS KT™**

The Delphis KT includes a laptop computer and all components housed on the modular and mobile KT cart. Adjustable shelves and workspace trays allow for customizable configurations. The Delphis KT has a wired connection between the laptop and UDS processor, but is Bluetooth™ capable as an option.



## **DELPHIS IP™**

The Delphis IP offers basic Urodynamics with software that is simple and easy to learn. An affordable, compact and portable unit mounted on an IV pole with a wireless connection to a laptop. This system is perfect for a clinic wanting to perform basic UDS studies and requiring a small footprint.




## **UROCAP™**

An innovative Uroflowmeter that provides accurate recordings which are easy to produce and interpret. The unit starts printing when flow is detected and stops when patient stops voiding. Automatic calculation of Uroflow Summaries. Available in a wired configuration with thermal printer or wireless to a computer.



## **NUWAV™ ULTRASOUND**

A system of ultrasound probes and unique software modules for bladder imaging, kidneys, pelvic floor or the prostate. Operates independently or integrated into a Laborie UDS device for a multi-function, multi-modality system.



## **UROSTYM™**

A biofeedback & stimulation unit that provides the most current protocols for the treatment of urinary incontinence. It generates reports quickly and easily that include complete patient information for analysis, treatment, and results. Urostrom can be added as an accessory to a UDS system or purchased as a stand-alone device.



## **A.R.M.**

Optimize your pelvic floor assessment with anorectal manometry. Our A.R.M. software provides the most accurate and comprehensive results for your diagnosis. A.R.M. can be added as an accessory to our UDS system or purchased as a stand-alone system.



U.S.A.

400 Avenue D, Suite 10  
Williston, VT 05495-7828  
Tel.: 1 800 522 6743  
1 802 857 1300  
Fax: 1 802 878 1122

INTERNATIONAL

6415 Northwest Drive  
Mississauga, ON L4V 1X1  
Tel.: +1 905 612 1170  
Fax: +1 905 612 9731

CANADA

2101 Boulevard Lapinière  
Brossard, Québec J4W 1L7  
Tél.: +1 450 671 5901  
Fax: +1 450 671 7182

EUROPE

Unit 5, Garonor Way  
Royal Portbury, Bristol  
United Kingdom BS20 7XE  
Tel: +44 1275 372332  
Fax: +44 1275 375145

[marketing@laborie.com](mailto:marketing@laborie.com) [www.laborie.com](http://www.laborie.com)